



10 IN 1 SUPER BUBBLE



- ★ Dig deep into pores and carry away dirt
- ★ Increases skin elasticity and brightness
- ★ Dispel acne and relieve skin sensitivity
- ★ Improve potholes and marks


Cosderma Aesthetica Pvt. Ltd

www.cosderma.org
www.cosdermamedicare.com


+91 9503 487 924 / +91 9987 686 006

Deep cleansing
skin system

Clean blackheads, acne,
worms; At the same time
add water, nutrients;
Narrow pores and
increase skin
transparency: Improve
skin.



Hydrogen and
Oxygen bubbles



Powerful
nine in one effect

Oxygen injection

EMS Nano
mesotherapy

Cold&Hot
skin handle

Ultrasonic
lead-in

BipolarHR

Skin analyser

Bubble Cleaning

Roller
speed lift

Oxygen
bubble

Ultrasonic
Probe

NINE POWERFUL HANDLES



Bubble Cleaning

A new generation of instruments, breaking through the technical level of traditional small bubble instruments, can open pores, soften keratin, make skin cleaner and deeper. Improve greatly the dark skin and blackheads, and make skin clean, white and tender.



Roller speed lift

The anion-cation EMS microcurrent is exchanged by positive and negative so that the cation is introduced into anion.



Ultrasonic lead-in

Ultrasonic vibrations of 1 to 3 million per minute frequency, open the pores, make the ingredients fast. Rapid penetration into the basal layer of the skin, allowing absorption rate more than 90%.



Oxygen injection

High pressure oxygen injection, particle penetration deep into the muscle. The skin layer enhances water retention by narrowing pores more. Easy to absorb.

Oxygen bubble

The formation of oxygen supply environment activated cells to enhance absorption, Leading the products into deep skin layer. Accompanied with the fruit acid products to achieve the result.



Cold&Hot skin handle

Fast skin qualitative, achieve the effect of anti-allergy and skin shaping



Ultrasonic lead-in

The skin scrubber uses high frequency vibration waves to soften skin. Keratinocytes deeply cleanse skin dirt.



Bipolar HR

The probe sends out tremors of up to six million per second. Radio waves, high frequency vibrations can make skin. The molecules in the skin generate heat from friction, which creates heat. The skin firming, improve the effect of wrinkles.



Wedinard Skin Diagnosis

Compatible with a smart Android system. Six major high-definition images analyze facial problems.


- Skin tone analysis
- Water and oil balance
- Fine line analysis
- Inflammation analysis
- Swatch analysis
- Pore analysis






Nanocrine uses nanoscale contacts, which can be as thin as more than 250 nanometers, and as short as 0.1 mm, just open the cuticle of the skin, without touching the dermis, which not only solves the problem of skin absorption, but also makes it painless and non-invasive!

Natural spread absorption




Natural absorption only penetrates into the skin layer of skin care molecules are difficult to penetrate and absorption is very low

Nanocrystalline nanocrystalline chips



The use of nanocrystals only through the epidermal barrier layer will not touch the blood vessels and nerves painless and non-invasive, and the absorption amount is greatly improved

EMS function



Microelectronic thermal conductivity technology can increase the absorption rate of the product by more than 40%

What changes it will bring to your skin?




Application: Quadrupole vacuum negative pressure radio frequency (high frequency) acts on the face.

The suction mode is recommended to adjust 1, the suction strength is recommended to adjust 1, and then manually adjust the suction time, it is recommended to adjust 04 or 06, and the operation is adjusted according to the customer's reaction.

Apply to the face and lift upwards according to the skin texture for 10 minutes on one side.

Specifications



- Voltage: 100v-120v 60Hz 200v-240v 50Hz
- Power: 200W
- HR frequency: 0.2-0.5 Mhz
- cold and hot scope :Cold 0-6C Hot 36-46C
- Ultrasound:1.1MHz
- Probe: 3 soft heads,8 hard heads
- Clean water mill:Positive pressure 4-5KgNegative pressure 96Kpa.Flows50L
- Display: 10.1inches Resolution:1024X600
- Solution bottle: A.Deep layer cleaning
 - B.Whitening and moisturizing
 - C.Deep cleaning
 - D.Instrument piping cleaning
- Skin scrubber vibration frequency: 26-30MHz
- Ems Beauty plastic head size:4X4mmX250Um
- RollerSpeedup:frequency200Hzmicrocurrent50-800(Ua)
- Displayscreen:10.1inches

Product photos





THANKS

Welcome to consult more!

www.cosderma.org / www.cosdermamedicare.com

+91 9503 487 924 / +91 9987 686 006

