

RENEULIFT

HIFU

Face Lifting | Wrinkle Removal | Body contouring

FMU BASED HIFU

SKIN TIGHTEN
&
SLIMMING

One Treatment Obvious Result



RENEULIFT

HIFU

Face Lifting | Wrinkle Removal | Body contouring

New Standard of MFU System

- No. 1 MFU proven by the patents
- MFU proven by cumulative sales and clinical cases
- Fastest treatment by Two-way round shooting and Auto Shot mode

What is MFU?

- | | | | |
|---|--|--|--|
| ■ Inner Skin Density Improvement Compared to HIFU | ■ Pores Count Improvement Compared to HIFU | ■ Skin Irregularity Improvement Compared to HIFU | ■ Jawline Improvement Compared to HIFU |
|---|--|--|--|

Faster shot speed

Efficacy improved by fastly repetitious shooting with 300 shots in 8 min.
Entertain more patients with faster treatment procedure.

Just 8 minutes Per 300 shot
MAXIMUM SPEED OF DOUBLO SERIES



RENEULIFT

HIFU

Face Lifting | Wrinkle Removal | Body contouring

Convenient treatment system

by Foot switch, by Hand piece




Foot switch



Hand piece

Intuitive GUI

Easy setting by intuitive button control




RENEULIFT


HIFU

Face Lifting | Wrinkle Removal | Body contouring

MFU treatment theory



Technology trend



RENEULIFT

HIFU

Face Lifting | Wrinkle Removal | Body contouring

Advantages of MFU over RF and laser

	MFH	RF	Laser
Number of Treatment	1 Time	4~5 Times	5 Times
Treatment Time	8 Min	30 Min	10~15 Min
Effect Durability	18 Months	1 Month	2~3Months
Depth	SMAS (3.0/4.5)	Below 3mm	1.5mm

For patients



RENEULIFT

HIFU

Face Lifting | Wrinkle Removal | Body contouring

Extremely Fast and safe **HIFU** solution
for face lifting & body contouring

Ultra-fast 3rd Generation
Transducers
300 shots in just 8 min.



- Face Lifting
- Neck treatment
- Eyebrow lifting
- Skin tightening and rejuvenation
- Body shaping
- Visible results at the first care.

No surgery No unavailability!
Painless effective and safe.



RENEULIFT


HIFU

Face Lifting | Wrinkle Removal | Body contouring




Product Description

HIFU (High intensity intensive ultrasound) energy is focused at the SMAS layer 3 to 4.5mm below the skin surface to produce a thermal coagulation zone smaller than 1 mm. At this zone, a high thermal reaction occurs at a temperature of 60°C or higher, which causes the SMAS layer contraction and collagen remodeling to skin tightening and lifting



It is unnecessary to worry about the skin injury because the power is across the epidermis. Moreover, it has the functions of quick pulling, counter compacting and quick wrinkle smoothing.




Application

1. Transducer M7 used to transmit ultrasound energy (high-intensity focused) to the epidermis tissue, the ultrasound to penetrate the skin to a depth of 1.5mm, is responsible for activating the skin epidermis in thinner tissue.

2. Transducer D4, used to release ultrasound energy into the dermis layer of 3.0mm, is responsible for activating the skin's dermal layer of collagen, effectively enhance the effect of the consolidation of the outline, but also improve large pores and reduce the appearance of wrinkles.

3. Transducer D7, used to transmit high-energy focus ultrasound to the skin subcutaneous tissue, the ultrasound to penetrate the skin to a depth of 4.5mm, direct subcutaneous SMAS layer, forming "heat coagulation" regional, targeted to thicker skin, such as cheeks, etc.



RENEULIFT

HIFU

Face Lifting | Wrinkle Removal | Body contouring



Product Effects

1. Skin Tightening and Rejuvenation
2. Non - Invasive Lifting for Face
3. Collagen Remodeling and SMAS Contraction



Specification

Type	Korea HIFU Ultrasound SMAS Lifting
Material	ABS
Power	110V to 240V 50ZH/60ZH
Frequency	7MHz 4MHz
Energy Power	0.1-2.0J
Mode	Pulse Digital output mode
Treatment cartridge	M7, D7, D4
Distance	Automatic adjustment
Length	5.0-2.5mm 5.0/step



RENEULIFT

HIFU

Face Lifting | Wrinkle Removal | Body contouring




HOW IT WORKS?

Reneulift HIFU (High Intensity Focused Ultrasound) is a non-surgical technology which offers brow lifting, jowl lifting, nasolabial fold reduction, periorbital wrinkle reduction and overall skin tightening and rejuvenation. The highly focused acoustic energy creates thermal coagulation zones at 3 different depths. Wound healing response results in the formation of new collagen thus providing longer term tightening of the skin. There is no down time associated with this non invasive procedure. HIFU has quickly become one of the most sought after non-surgical face lifting treatments with the need for only one session, and results which can last up to 18 months



INTERFACE



User friendly interface with "Professional" and "Intelligent" 2 working models



THANKS

Welcome to consult more!



www.cosedermapeel.com www.cosdermamedicare.com

+91 9503 487 924 / +91 9987 686 006